

Improving the productivity of your chemistry

Findenser[™] super air condenser

MIR

A HEAR

Patented

Replaces the need for water-cooled condensers in over 95% of all common chemistry applications

• Save water

Prevent flooding

www.radleys.com

Findenser - minimises risk of flooding in the laboratory

Findenser requires no running water to operate. Water is a precious resource. It makes little economic or environmental sense to waste thousands of litres just to cool a single condenser.

> Ground glass socket

 Accepts standard tapered cones

Choice of sizes

Easy clamping

• Use with standard laboratory clamps Clamping kits available

Anti-roll design

What are the benefits?

- No risk of flooding from running water
- Eliminates water purchase and disposal costs
- For solvent volumes from 5ml up to 1 litre
- Helps meet sustainable water reduction targets


Finned aluminium jacket

- Excellent thermal conductivity
- High performance air-cooling Chemical-resistant
 - anodised surface



The cost of running a	1 tap running at 2.5 litres per minute			
water condenser	4 hours	8 hours	24 hours	
Water consumption per day (litres)	600	1200	3600	
Cost per day	£1	£2	£6	
Cost per month, 20 working days	£20	£40	£120	
Cost per year	£240	£480	£1,440	

Based on a flow of 2.5 litres per minute, at a cost of 0.16p per litre for supply and disposal of water. Price from Veolia Water 2015 commercial pricing



2

Findenser - performs as well as a water condenser

Findenser replaces water-cooled condensers in over 95% of common chemistry applications

Findenser performance testing

A range of solvents, in identical flasks and set-ups, were tested with a Findenser, water condenser and air condenser to record solvent loss by weight.


Compared with an air condenser

For synthesis with low boiling point solvents. Findenser showed a significant improvement in solvent retention. With acetone or DCM the reaction boiled dry when using an air condenser, yet Findenser retained 95% of the solvent under the same conditions.

For synthesis with medium boiling point solvents, Findenser delivered improved solvent retention particularly with larger volumes and high temperatures.

Compared with a water condenser

Under identical conditions, Findenser retained solvent to the same level (within the limits of experimental error) as a water condenser with a range of solvents. Diethyl ether is not recommended for use with Findenser.



How does Findenser work?

- Findenser comprises an internal glass condenser and an external, finned aluminium jacket, between which a small amount of water is permanently sealed.
- The glass condenser design has a greater internal surface area than traditional condensers, increasing heat transfer capacity.
- The finned jacket fits around the glass condenser, further increasing the external surface area.
- The result is a 'SUPER air condenser'.


Choice of B14, B19, B24 and B29 joint sizes





Time to switch





Findenser Mini cross section

Findenser ordering information

Findenser RR31100 RR31102 RR31104	Findenser B24 Cone, B24 Socket - 400mm long Findenser B29 Cone, B24 Socket - 400mm long Findenser B19 Cone, B19 Socket - 400mm long				÷	► 400 H ► 400 H (0 /
Findenser Mir RR31105 RR31107 RR31109	i Findenser B24 Cone, B24 Socket - 275mm long Findenser B19 Cone, B19 Socket - 275mm long Findenser B14 Cone, B14 Socket - 275mm long		RR71510 Rod ar for Carousel		RR71516 Rod ar for Carouse	
Accessories RR31200 RR31202 RR31204 RR31206 RR31208 RR139139	PTFE Adapter B24 Socket to B29 Cone PTFE Adapter B19 Socket to B24 Cone PTFE Adapter B19 Socket to B29 Cone PTFE Adapter B14 Socket to B19 Cone PTFE Adapter B14 Socket to B24 Cone Right Angle Adapter Ordinary B24 + GL14 + fitting	15			0	
RR139140 RR139141 RR139142	Right Angle Adapter Ordinary B19 + GL14 + fitting Right Angle Adapter Ordinary B19 + GL14 + fitting Enlarging Adapter Ordinary B24 to B29	js	RR71125 Support Rod S RR71127 Support Rod S		RR71120 Support Rod Hotplate Adapter	
Clamps and A RR31210 RR31212 RR31214 RR31216 RR31218	ccessories Retort Stand, 2 Position Base + 12x750mm SS R Retort Stand, 2 Position Base Support Rod SS 12x600mm No Thread Support Rod SS 12x750mm No Thread Support Rod SS 12x1000mm No Thread	od			34 140	10 m
RR71110 RR71115 RR71120	Refort Clamp to 85mm Boss Head to 16mm Support Rod Hotplate Adapter (extension plate)		RR71512 Stand, Ro	d and Clamp Kit	RR71514 Stand, Ro	d and 2 x Clamp Kit
RR71125 RR71127 Rod and Clan	Support Rod SS 13x340mm M10 Support Rod SS 13x500mm M10		1		1	
RR71510	Rod and Clamp Kit for Carousel Hotplate 1 x RR71120 - Support Rod Hotplate Adapter 1 x RR71127 - Support Rod SS 13x500mm M10 1 x RR71110 - Retort Clamp to 85mm 1 x RR71115 - Boss Head to 16mm		RR31210 Retort Stand		RR31212 Retort Stand	4
RR71516	Rod and 2 x Clamp Kit for Carousel Hotplate 1 x RR71120 - Support Rod Hotplate Adapter 1 x RR71127 - Support Rod SS 13x500mm M10 2 x RR71110 - Retort Clamp to 85mm 2 x RR71115 - Boss Head to 16mm		+ 12x750mm		2 Position Base	
RR71512	Stand, Rod and Clamp Kit 1 x RR31210 - Retort Stand, 2 Position Base and 12x750mm SS Rod 1 x RR71110 - Retort Clamp to 85mm 1 x RR71115 - Boss Head to 16mm	**	RR31216 Suppor	rt Rod SS 12x600mm I rt Rod SS 12x750mm I t Rod SS 12x1000mm	No Thread	
RR71514	Stand, Rod and 2 x Clamp Kit 1 x RR31210 - Retort Stand, 2 Position Base and 12x750mm SS Rod 2 x RR71110 - Retort Clamp to 85mm 2 x RR71115 - Boss Head to 16mm		S.	- Carlos	1008	T
	at both Findenser and your flask are securely clamped.		RR71110 Retort	RR71115 Boss	RR139139/140	RR31200 PTFE
Technical Length x Diam	Specifications eter 275mm x 72mm / 400mm x 72mm		Clamp to 85mm	Head to 16mm	Right Angle Adapters	Adapter B24 Socket to B29 Cone
Joints Weight Inner or wetter	B14, B19, B24 & B29 670g / 1200g					

www.radleys.com

Finned jacket

Plug and Seal

Max solvent bpt

Operating temp. range

Thermofluid

Oven drying

Dishwasher

Autoclaving

World leaders in innovative productivity tools for chemists

Anodised aluminium

Max temp. 60°C Max temp. 50°C

Do not autoclave 0°C to 60°C

Water

155°C

Acetal and HT silicone

rodeys innovations for chemistry

> Shire Hill, Saffron Walden, Essex, CB11 3AZ. United Kingdom. t: +44 1799 513320 f: +44 1799 513283 e: sales@radleys.co.uk