pipetman®

4

THE PIPETTING STANDARD RECOGNIZED RELIABILITY AND PROVEN DURABILITY COMBINED WITH EASY PIPETTING AND HANDLING



Legendary Performance

Accurate and precise, PIPETMAN[®] air-displacement pipettes are recognized as one of the world standard. They come in a large range of volumes, including single and multichannel models, allowing you to handle various sample types with confidence.

Easy Pipetting

PIPETMAN offers low pipetting forces while preserving standard pipetting habits. PIPETMAN allows hours of comfortable and easy pipetting.

Built-to-Last and Convenient

Made of PVDF and stainless steel, PIPETMAN offers premium durability. Built-to-last, it is the economical solution for years of pipetting.

Simple to Maintain

Routine cleaning and maintenance are all that is required to keep PIPETMAN in top condition for years. The PIPETMAN design allows for easy access to userserviceable parts, along with a fully autoclavable tip holder and ejector.





Choose Your Pipette Stand

Different pipette stands and holders allow you to store pipettes vertically, protecting them from contamination.



Pipette Stands				
Description	Part Number			
CARROUSEL™ Pipette Stand (7 pipettes)	F161401			
TRIO™ Pipette Stand (3 pipettes)	F161405			
Universal Multichannel Stand	F161417			
SINGLE™ Pipette Holder	F161406			

Choose Your Pipette Kits



PIPETMAN single channel models are available in kits, offering the essential pipetting tools to start an experiment right away.

PIPETMAN* Kits					
Model	Part Number	Models			
PIPETMAN® 4-Pipette Kit	F167360	P2 P20 P200 P1000			
PIPETMAN® Starter Kit	F167900	P20 P200 P1000			
PIPETMAN® Micro-volume Kit	F167800	P2 P10 P100			
PIPETMAN® Large-volume Kit	F167920	P5000 P10mL			

PIPETMAN [®] Single Channel					
Model	PIPETMAN [®] DIAMOND Tips		Part Number	Volume	
P2	D10 DL10	DF10ST DFL10ST	F144054M	0.2-2 μL	
P10	D10 DL10	DF10ST DFL10ST	F144055M	1-10 µL	
P20	D200	DF30ST	F144056M	2-20 μL	
P100	D200	DF100ST	F144057M	10-100 μL	
P200	D200 D300	DF200ST DF300ST	F144058M	20-200 μL	
P1000	D1000 D1200	DF1000ST DF1200ST	F144059M	100-1000 μL	
P5000*	D5000		F144066	500-5000 μL	
P10mL*	D10mL		F144067	1-10 mL	
PIPETMAN* Multichannel					
P8x10 P12x10	D10 DL10	DF10ST DFL10ST	F144068 F144069	1-10 μL	
P8x20 P12x20	DL10 D200	DFL10ST DF30ST	F144070 F144071	2-20 μL	
P8x200 P12x200	D200 D300	DF200ST DF300ST	F144072 F144073	20-200 μL	
P8x300 P12x300	D200 D300	DF200ST DF300ST	F144074 F144075	30-300 μL	

Models equipped with a tip ejector are only sold with a stainless steel tip ejector. Plastic tip ejector is available as an accessory. Refer to *PIPETMAN* User's Guide LT801122*. *P5000 and P10mL do not have tip ejectors.

Gilson maximum permissible errors are guaranteed only when PIPETMAN $^{\rm e}$ pipettes are used with the recommended PIPETMAN DIAMOND Tips.



PIPETMAN DIAMOND Tips are designed to fit your PIPETMAN perfectly and guarantee maximum accuracy and precision.

For ordering information, visit www.gilson.com

gilson.com/contactus

LT800591/E [@2023 Gilson S.A.S. All rights reserved. All trademarks are the property of Gilson, Inc. and its subsidiaries.

