Agilia Rolling Stands for Increased Mobility

Multi-Channel Rolling Stand

Specifically designed for the mobility of Link Agilia ensuring a safe centre of gravity.



Doubled pole to ensure well positioned centre of gravity

5 doubled castors with brakes

Robust weighted central base

Tube diameter: 34 mm

Telescopic bag hanger (170 to 220 cm) with 6 hooks for up to 6 kg in total

Can accommodate 8 Agilia pumps and Link Agilia racks up to a total of 21.8 kg

Art. Nº: Z073160





Fresenius Kabi AG 61346 Bad Homburg Germany Phone +49 (0)6172 686-0 www.fresenius-kabi.com

🕨 Link Agilia range

Offers a robust column to improve the organisation of 4, 6 or 8 Agilia pumps



Folding Handle

Optional accessory to facilitate handling and transport. Art. N°: 2073117





Optional accessory to support light objects, 2 kg max. Art. N°: Z073113



U-support

Optional accessory to fix light IV devices, 2 kg max. **Art. N°: Z073116**



Agilia is a registered trademark by Fresenius Kabi in selected countries.

Agilia Rolling Stands for Increased Mobility

Infusion Pump Rolling Stand

Designed for the safe and stable mobility of 1, 2 or 3 infusion pumps.



All-in-one product 5 doubled castors with brakes Robust weighted central base Tube diameter: 25 mm Telescopic bag hanger (170 to 220 cm) with 4 hooks for up to 2 kg each Integrated handle Can accommodate up to 3 Agilia pumps (or up to 6.8 kg) Integrated handle

Art. Nº: Z073150

Due to our policy of continuous product development as well as changes in standards, the features described are subject to change. Manufactured by Provita. Distributed by Fresenius Vial.



Fresenius Kabi AG 61346 Bad Homburg Germany Phone +49 (0)6172 686-0 www.fresenius-kabi.com

From General Wards to Intensive Care units,

the Agilia family offers a complete range of infusion pumps, disposables and accessories to cover all major IV drug administration needs.



Agilia is a registered trademark by Fresenius Kabi in selected countries.

Agilia Rolling Stand for Increased Mobility

Twin Link Rolling Stand

Specifically designed for the mobility of 2 Link.



Robust weighted central base to ensure stability 4 doubled antistatic castors with brakes Tubes diameters: 38 mm Device weight: 46 kg Telescopic bag hanger (184 to 220 cm) with 6 hooks 2 additional hooks on the Twin Link structure Up to 2 Link Agilia and 16 Agilia pumps Max devices weight on each poles: 25 kg (50 kg total) Max container weight on each hook: 2 kg (16 kg total)

Art. Nº: Z073170

Due to our policy of continuous product development as well as changes in standards, the features described are subject to change. Manufactured by Provita. Distributed by Fresenius Vial.



Fresenius Kabi AG 61346 Bad Homburg Germany Phone +49 (0)6172 686-0 www.fresenius-kabi.com



🕨 Link Agilia range

Offers a robust column to improve the organisation of 4, 6 or 8 Agilia pump



Agilia is a registered trademark by Fresenius Kabi in selected countries.