









PIPETTE TIP COMPATIBILITY CHART

*Pipette tip series marked with asterisk are not compatible with multichannel pipettors.

Compatibility is determined by fit of tip onto the pipettor, and the ability to meet manufacturer's specification when pipettors are calibrated with tip.

Pipette Tip Compatibility Chart

Pipettor Brand/Model	Biohit mLine™	Biohit Proline™	Biohit Proline Plus™	Brand Transferpette S™	Brand Transferpette Electronic™	Capp™	Eppendorf® Reference™	Eppendorf® Research™	Eppendorf® Research Plus™
2040 Series – 10 µl Micro Tip 		0.5-10 µl		0.5-10 µl		0.5-10 µl	0.1-2.5 µl 0.5-10 µl	0.1-2.5 µl 0.5-10 µl	0.1-2.5 µl 0.5-10 µl
2340 Series – 10 µl Extended Length Tip 	M10	0.5-10 µl		0.5-10 µl		0.5-10 µl	0.1-2.5 µl 0.5-10 µl 2-20µl	0.1-2.5 µl 0.5-10 µl	0.1-2.5 µl 0.5-10 µl
2140 Series – 10 µl E-Style Tip 	M10	0.5-10 µl		0.5-10 µl		0.5-10 µl	0.1-2.5 µl 0.5-10 µl 2-20µl	0.1-2.5 µl 0.5-10 µl	0.1-2.5 µl 0.5-10 µl
2100 Series – 200 µl Universal Tip 	M100 M200		10-100 µl 20-200 µl	2-20 µl 10-100 µl 20-200 µl	15-300 µl	5-50 µl 10-100 µl 30-300 µl	2-20 10-100 50-200 µl	2-20 µl 10-100 µl 20-200 µl 30-300 µl	2-20 µl 10-100 µl
2016 Series – 200 µl Extended Length Gel Tip* 	M100 M200		10-100 µl 20-200 µl	2-20 µl 10-100 µl 20-200 µl	15-300 µl	5-50 µl 10-100 µl 30-300 µl	2-20 10-100 50-200 µl	2-20 µl 10-100 µl 20-200 µl 30-300 µl	2-20 µl 10-100 µl
2150 Series – 200 µl Extended Length Tip* 	M100 M200		10-100 µl 20-200 µl	2-20 µl 10-100 µl 20-200 µl	15-300 µl	10-100 µl 30-300 µl	2-20 10-100 50-200 µl	2-20 µl 10-100 µl 20-200 µl 30-300 µl	2-20 µl 10-100 µl
2090 Series – 300 µl Universal Tip 	M100 M200 M300		10-100 µl 20-200 µl 30-300µl	20-200 µl	15-300 µl	5-50 µl 10-100 µl 30-300 µl	50-200 µl	20-200 µl 30-300 µl	50-200 µl
2110 Series – 1000 µl Traditional Shaped Tip* 	M1000	200-1000 µl		100-1000 µl		100-1000 µl	100-1000 µl	100-1000 µl	100-1000 µl
2160 Series – 1000 µl Universal Tip* 	M1000	200-1000 µl		100-1000 µl		100-1000 µl	100-1000 µl	100-1000 µl	100-1000 µl
2370 Series – 1000-1250 µl Extended Length Tip 	M1000	200-1000 µl		100-1000 µl		100-1000 µl	100-1000 µl	100-1000 µl	100-1000 µl
2250 Series – 5000 µl Tip 									

Eppendorf® Xplorer Plus™	Finnpipette™	Finnpipette F2™	Finnpipette Digital™	Finnpipette F1™	Gilson Pipetman™	Gilson Pipetman Ultra™	Hamilton™	Labnet BioPette E™	Nichiryo Nichipet EX™	Nichiryo Oxford Benchmark™	Nichiryo Oxford Multimate™	Socorex Calibra 822™	VWR Ergonomic High Performance™	VWR Ultra High Performance™
	0.5-10 µl 5-50 µl		1-10 µl 0.5-10 µl	0.2-2 µl 1-10 µl	P2 P10		0.2-2 µl 1-10 µl		0.5-10 µl	0.1-2 µl	0.5-10 µl	1-10 µl		0.1-2 µl 0.5-10 µl
	0.5-10 µl 5-50 µl		1-10 µl 0.5-10 µl	0.2-2 µl 1-10 µl	P2 P10		0.2-2 µl 1-10 µl		0.5-10 µl	0.1-2 µl	0.5-10 µl	1-10 µl		0.1-2 µl 0.5-10 µl
	0.5-10 µl 5-50 µl		1-10 µl	0.2-2 µl 1-10 µl	P2 P10		0.2-2 µl 1-10 µl		0.5-10 µl	0.1-2 µl	0.5-10 µl	1-10 µl		0.1-2 µl 0.5-10 µl
	2-20 µl 5-50 µl 20-200 µl 30-300 µl	2 - 20µl 5 - 50µl 10 - 100µl 20 - 200µl	5-50 µl 20-200 µl 30-300 µl	2-20 µl 5-50 µl 20-200 µl	P20 P100 P200	U20 U200	2.5-25 µl 10-100 µl 30-300 µl		2-20 µl 10-100 µl 20-200 µl	2-20 µl	5-50 µl 50-300 µl	10-100 µl 20-200 µl	2-20 µl 20-200 µl	2-20 µl 10-100 µl 20-200 µl
	2-20 µl 5-50 µl 20-200 µl 30-300 µl		30-300 µl		P20 P100 P200	U20 U200	2.5-25 µl 10-100 µl 30-300 µl		2-20 µl 10-100 µl 20-200 µl	2-20 µl		10-100 µl 20-200 µl	2-20 µl 20-200 µl	2-20 µl 10-100 µl 20-200 µl
	5-50 µl 20-200 µl 30-300 µl		30-300 µl		P100 P200	U200	2.5-25 µl 10-100 µl 30-300 µl		10-100 µl 20-200 µl			10-100 µl 20-200 µl	20-200 µl	10-100 µl 20-200 µl
15- 300µl	5-50 µl 20-200 µl 30-300 µl		30-300 µl 50-300µl	30-300 µl	P200	U200	2.5-25 µl 10-100 µl 30-300 µl		20-200 µl		50-300 µl	20-200 µl	2-20 µl 20-200 µl	2-20 µl 10-100 µl 20-200 µl
50- 1000µl	100- 1000 µl 200- 1000 µl	100 - 1000µl	100-1200 µl	100 - 1000µl	P1000		100-1000 µl	100- 1200 µl	100-1000 µl			100-1000 µl	100- 1000 µl	100-1000 µl
50- 1000µl	100- 1000 µl 200- 1000 µl	100 - 1000µl	100-1200 µl	100 - 1000µl	P1000		100-1000 µl	100- 1200 µl	100-1000 µl			100-1000 µl	100- 1000 µl	100-1000 µl
50- 1000µl	100- 1000 µl 200- 1000 µl	100 - 1000µl	100-1200 µl	100 - 1000µl	P1000		100-1000 µl	100- 1200 µl	100-1000 µl			100-1000 µl	100- 1000 µl	100-1000 µl
					P5000	U5000								

Pipette Tip Compatibility Chart

