

PIPETTE CALIBRATION 2026



Accredited calibration, *As Found* and *As Returned*:

Ramcon is accredited according to ISO17025:2018 and performs calibration of your pipettes. It consists of 3x10 measurements, for both *As Found* and *As Returned*. The calibration is performed at 100%, 50% and 10% of the pipettes nominal volume.

A calibration is performed when we first receive the pipette, *As Found*. Repairs and adjustments are carried out as needed. After adjustment, a new calibration is made, *As Returned*. You will receive a calibration report with the results from both *As Found* and *As Returned*.

Calibration *As Found* and *As Returned*, or only *As Returned*:

2x5 measurements of maximum volume and minimum volume, both for *As Found* and *As Returned*. A calibration is performed when we first receive the pipette, *As Found*. Any repair and adjustment are carried out. After adjustment, a new calibration is made, *As Returned*. You will receive a calibration report with the results from *As Found* and one from *As Returned*.

All our internal methods are based on ISO 8655-2:2022.

Order calibration: service@ramcon.se



Typ av service Type of service	Basic	Combo	Advanced	Combo ACC	Advanced ACC
According to Internal Method ISO 17025	x	x	x	✓	✓
As found calibration	x	✓	✓	✓	✓
As returned calibration	✓	✓	✓	✓	✓
Cleaning	✓	✓	✓	✓	✓
Maintenance service	✓	✓	✓	✓	✓
Replacement of spare parts at no extra cost	x	x	✓	x	✓
Discount on parts used	x	10 %	100 %	10 %	100 %
Free return shipping	✓	✓	✓	✓	✓

